



## Cohesive Bandages Tan 2.5cm

**\$7.19**

In stock  
Only **997** left  
**SKU** cohesive-  
tan-2.5



 Add to Wish List

 Add to Compare

- Great for the kitchen
- Use as a strapping tape for sport
- Reusable



### Details

## Cohesive Bandage Tan 2.5cm

*Maintains compression and is great for the workplace or for physical exercise*

- Great for the kitchen
- Use as a strapping tape for sport
- Tan coloured to be more discrete
- Reusable

These cohesive bandages are a great addition to your first aid kit as they are one of your most versatile bandages.

Cohesive bandages require no adhesive, clips or pins, but will hold to itself and maintain great compression.



With hand-tear technology embedded into the bandage, no additional cutting materials like scissors are needed.

This bandage is so lightweight, breathable and comfortable that it is ideal for your workplace.

These non-sterile bandages can act as a secondary dressing.

Cohesive bandages with their tension retention can offer support for weak or injured joints e.g. sprained ankle, knee ligament injury

Unlike other strapping tapes, these cohesive bandages maintain their tension meaning that you do not have to make any adjustments when you apply

- Applying the bandage with stretch gives extra tension
- Applying the bandage without stretch provides light support
- Extra layers can of bandage can provide more support

Dimensions: 2.5cm x 4.5m

**Contents:**

1 x roll