



Fevermates - Stick On Thermometers Temperature Indicator (8 Pieces)

\$7.69

Out of stock
SKU FMSTICK01

 Add to Wish List

 Add to Compare

This innovative fever indicator allows you to visually monitor your children's temperature up to 48 hours with NO FUSS.

Fast, safe and effective Kids love to wear them

Hypo-allergenic

Latex-Free



Details

Fuss Free Fever Monitoring for 48 hours

Sick kids and thermometers don't mix - so use a Fevermate for up to 48 hours to monitor what's happening.

The convenient, fuss free stick-on fever and temperature indicators that allow you to continuously monitor your children's fever or temperature for up to 48 hours at a glance.

This innovative fever and temperature indicator allows you to visually monitor your children's temperature up to 48 hours with NO FUSS.

- Fast, safe and effective
- Kids love to wear them



- Hypo-allergenic
- Latex-Free

Makes reading temperature easy and fuss-free with easy-to-read clear numbers that display up when your child's temperature is reached. Reading thermometers is hard and kids often don't like the thermometer in their mouth or under their arm-pit. Works with adults too!

Contents: 8 Pack of Fevermates stick-on animal temperature indicators & Instructions

Complies with Health & Safety requirements: CE 0197.

ARTG Identifier: 232932. Included on the ARTG as a medical device. (ARTG = Australian Register of Therapeutic Goods).

Warning! Item is medical device and not a toy. Possible choking hazard! If applying to a child 36 months or younger, must be supervised at all times and not left unattended with item. Always use a thermometer to confirm temperature reading. If you have concerns, seek medical advice. Keep packaging and contents away from children.