




## Aerated Eye Wash Unit Free Standing, Pedestal MTD

**\$836.41**

In stock  
Only **1000** left  
**SKU** 99692

 [Add to Wish List](#)

 [Add to Compare](#)





## Details

Speakman SE-545 features a unique panic bar activator on top of the unit, making it easy to find and allows for either right or left hand activation. The soft aerated sprays allow users to comfortably flush their eyes for the ANSI recommended 15 minutes.

The SE-545 is pedestal mounted unit and comes complete with a stainless steel rectangular bowl, a trigger bar activator and is compliant to the ANSI Z358.1 standard.

## Features:

- Pedestal mounted
- Stainless steel rectangular bowl
- Aerated spray outlets
- Trigger bar activating system
- Meets ANSI Z358.1

## Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

## Related Products

Check items to add to the cart or select all