



Blood Pressure Monitor DuraShock DS66 Aneroid.

\$289.08

In stock
Only **1000** left
SKU 99650

 [Add to Wish List](#)

 [Add to Compare](#)

Details

DuraShock DS66 maintains accurate blood pressure readings with lower lifetime cost with gear-free DuraShock technology. The DS66 Trigger features an innovative trigger-release design for easy operation.

Features:

- Shock resistant to maintain accuracy with lower lifetime cost. it can withstand up to 30" drop and remain in calibration.
- Soft comfort grip.
- Jewel movement for long life.
- Laser-engraved dial face for accuracy.
- Lighter weight than traditional models.
- Premium aesthetics.
- Certified accuracy to +/- 3mmHg.
- Trigger deflation valve for smooth, one-finger operation.
- Latex free for safety.
- Durable one-piece Adult cuff.
- Polyester zip case.



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all