



## Blood Pressure Monitor Erka Kobold Professional One handed EK223

**\$243.31**

In stock  
Only **1000** left  
**SKU** 400114

 [Add to Wish List](#)

 [Add to Compare](#)



### Details

The ERKA Kobold is a robust and reliable partner for all professional users.

Thanks to the ergonomic, balanced design combined with optimisation for low weight and the highly sensitive threaded discharge valve, maximum ease of use is guaranteed in the daily routine of the medical professional – for every measurement.

Due to its ergonomic design and the sensitive screw-on bleeder valve, the ERKA Kobold is very easy to use. The robust case is made of anodised aluminium for optimised surface hygiene.

The 48mm scale makes the ERKA Kobold easy to read. The manual blood pressure monitor is equipped with the ERKA snap closure which makes it possible to change the cuff within a few seconds.



### Contains:

- Manual blood pressure monitor ERKA Kobold
- Connecting hose with ERKA snap closure
- Upper arm cuff in size 4 (arm circumference of 27 - 35 cm)
- Case

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all