



Blood Pressure Monitor Set Liberty Adult with Aneroid Spyhg & Stethoscope

\$108.25

In stock
Only **1000** left
SKU 400602

 [Add to Wish List](#)

 [Add to Compare](#)

Details

- Black enamel 300mm Hg no stop pin manometer
- D-Bar self adjusting velcro cuff for one handed application.
- Complete with bladder, bulb and valve.
- Chrome plated inflation valve.
- Stethoscope attached.
- Zipped leatherette case.
- Complete with instructions.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all