



Conforming Gauze Bandage Light Support

\$2.03

In stock
Only **999** left
SKU WEB-CGB

Size

-- Please Select -- ▼

♥ Add to Wish List

▬ Add to Compare

Details

Conforming Gauze Bandages are constructed from a blend of polyester yarn and rayon. These bandages are soft, light weight and porous, giving great comfort and high levels of air circulation.

Options:

- 2.5cm x 1.8m
- 10cm x 1.8m
- 15cm x 1.8m

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.