



## HeartSine Samaritan Pad 350P Defibrillator & Carry Case

**\$1,986.17**

In stock  
**SKU** 400687

 Add to Wish List

 Add to Compare

### Details

Sudden Cardiac Arrest (SCA) often strikes without warning. There's little time to react and no time to think, so it's important to have the HeartSine Samaritan PAD 350P defibrillator on hand.

The 350P is a clinically advanced, lightweight, compact, easy-to-use, and portable AED. Audio and visual prompts guide the user step by step through the rescue process. From trained first aiders to professional responders, the 350P is designed so that anyone, anywhere, anytime can save a life.

- Compact and lightweight portable AED
- Easy to use
- HeartSine's unique biphasic waveform.

The Samaritan PAD has a rating of IP56, providing the 350P defibrillator with the highest level of protection from dust or moisture ingress in the industry.

The 350P uses HeartSine's practical Pad-Pak, which incorporates the battery and electrodes in a single-use cartridge – meaning only one expiration date to monitor and only one item to replace after use.

The 350P is suitable for use on adults and children. (For children younger than 8 years or weighing less than 25kgs a Pediatric-Pak is available).



---

## Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.