




HeartSine Samaritan Pad 500P Defibrillator & Carry Case

\$2,643.81

In stock
SKU 400474

 [Add to Wish List](#)

 [Add to Compare](#)



Details

- ICG-based CPR Advisor gives the rescuer feedback on compressions based on actual cardiac output: "Push Harder," "Push Faster," "Push Slower" or reassurance that they are performing "Good Compressions."
- Compact and lightweight, 28–56% smaller than other portable AEDs on the market and weighing only 1.1kgs.
- Easy to use. Audio and visual prompts guide the user step by step through the rescue process. Three large, easy-to-see icons illuminate: "Attach pads," "Do not touch the patient," and "It is safe to touch the patient" to facilitate the rescue process.
- Ophisticated SCOPE biphasic waveform technology adjusts patient impedance to maximize effective shock delivery.
- HeartSine's practical Pad-Pak combines battery and electrodes into one cartridge – one date to monitor and only one item to replace on expiration or after use.
- Automatic self-check.
- The 500P defibrillator has a rating of IP56, providing the highest level of protection from dust and moisture ingress.
- Suitable for adults and children.
- Controls and electrodes are accessible even if the HeartSine Samaritan PAD is still in the case.
- Integrated metronome prompts the CPR chest compression rate.
- Events are recorded, with the time, date, ECG trace and duration, shock delivery information and CPR intervals.



- CE Marked.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.