



Hospital Crepe Bandage 5cm

\$2.76

In stock
Only **1000** left
SKU 4111

 [Add to Wish List](#)

 [Add to Compare](#)



Details

Handy Crepe is a high quality white cotton crepe bandage. Clips included.

Features:

- Non-fraying edges
- Clip to secure bandage included
- Cool and comfortable bandaging



Indications:

- Dressing retention where light compression is required
- Post-operative use

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.