



Liberty Mini Fibre Optic Otoscope (Auroscope)

\$155.18

In stock
Only **998** left
SKU 400122

 [Add to Wish List](#)

 [Add to Compare](#)



Details

- New improved switch and bayonet fitting for more secure attachment of specula.
- Fibre Optic.
- Bright 2.5 Volt F.O. Light source.
- Battery handle (2 x AA size batteries).
- Plastic tube with reusable ear specula.
- Also fits Heine disposable tips.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all