



## Liberty Pocket Otoscope & Ophthalmoscope Diagnostic Set

**\$384.79**

In stock  
Only **1000** left  
**SKU** 400605

 Add to Wish List

 Add to Compare



### Details

The Liberty Pocket Otoscope Ophthalmoscope Diagnostic Set is super convenient, portable and value-packed with strong metal pocket clips. An ophthalmoscope giving 6 diagnostic beams, and a closed head otoscope, complete with 4 extra specula, use both in the consulting room and when on the move. Supplied in a hard protective carry case.

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other



qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

**Related Products**

Check items to add to the cart or select all