



## Little Anne QCPR

**\$667.81**

In stock  
Only **1000** left  
**SKU** 123-01050

 Add to Wish List

 Add to Compare



### Details

## Laerdal CPR-D (Defibrillator) Manikin AED Little Anne

### Life-like anatomy

Perform nose-pinch, head tilt, chin lift and jaw thrust to see chest rise. Anatomically correct oral and nasal passages.

### Correct compression force

See and hear that compressions are done correctly. A built-in clicker signals when the correct compression depth is achieved.

### Quick & easy maintenance

Clean-up is easy with disposable lungs and removable faces. Made with Laerdal's patented design and materials; it's built to last.

Little Anne QCPR allows an instructor to monitor up to six learners at the same time (using our QCPR Training app) or up to 42 learners (with QCPR Classroom) from the palm of their hand. The Little Anne QCPR Training and QCPR Classroom Apps provide an easy to view performance summary that highlight which students may require further guidance.

At a glance you can pick up who is not compressing deeply enough, who is not compressing at the correct rate and who is overventilating.



This product includes:

Manikin

1 additional manikin face

1 additional airway

6 manikin wipes

Jacket for manikin

Soft pack carrying case/training mat

2 AA Batteries

1 O-Ring

Little Anne Manikin from Laerdal designed for use with the Laerdal AED2 and AED3 training defibrillators.

Training defibrillators sold separately - this listing is for the manikin only.