



Med Sled Evacuation Stretcher with 3 Cross Straps & Tether

\$1,324.12

In stock
Only **999** left
SKU 280105

 [Add to Wish List](#)

 [Add to Compare](#)



Details

When timing is critical, Med Sled Evacuation enables employees to safely transport and evacuate patients from a single or multilevel building.

The stairwell braking system and integral patient support provides control of descent with patients two to three times their own weight. No lifting is required when using the roll and drag method.

Each Sled has eight handles for lifting, three securing cross straps and a tether and stairwell braking system. The Med Sled Evacuation rolls up and fits into the supplied compact bag with shoulder strap.

Toddler and Infant padded inserts to fit within the MS36 are available. Wall unit storage devices available for sleds and inserts.

- Dimensions: 91 x 221cm
- Weight: 2.7kg

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

[Check items to add to the cart or select all](#)

